Vet Expert Rapid Test FPV/FCoV/Giardia Ag

EN PRINCIPLES

The Vet Expert FPV/FCoV/Giardia Ag is an immunochromatographic test for the qualitative detection of Feline Parvovirus (FPV), Feline Coronavirus (FCoV) and *Giardia duodenalis* antigens in feline feces. The test is performed on a test plate indicating test line (T) and control line (C). The control line (C) is used for procedural control and should always appear if the test procedure is being performed properly and the test reagents of the control line are working. A colored test line (T) will be visible in the result window if the target antigens are present in the specimen.

Sensitivity: FPV 97.1%/FCoV 92.9%/Giardia 99.9%
Specificity: FPV 98.5%/FCoV 98.1%/Giardia 99.9%

MATERIALS PROVIDED (2/5/10 TESTS/KIT)

- 2/5/10 Vet Expert FPV/FCoV/Giardia Ag test plates
- 2/5/10 tubes with diluent for FPV/FCoV (P)
- 2/5/10 tubes with diluent for Giardia (G)
- 4/10/20 collection swabs
- 4/10/20 disposable pipettes
- Instruction for use

PRECAUTIONS:

- For veterinary diagnostics use only.
- Perform the test immediately after opening the pouch.
- Strict adherence to the instruction is necessary.
- The test components are single-use only.
- Do not touch the membrane in the result window.
- Do not use the test kit beyond the stated expiration date marked on the package label.
- Do not use the test kit if the pouch is damaged or the box seal is broken.
- Do not mix components from different lot numbers because the components in this kit have been quality control tested as a standard batch unit.
- Handle with the sample as potentially infectious.
- Decontaminate and dispose of all samples, used kits and potentially contaminated materials safely in accordance with national and local regulatory requirements.
- Consider the test results as invalid after the specified read out time.

STORAGE OF THE TEST KIT

The test kit can be stored between 2-30°C. The test is stable within the expiration date marked on the packaging. DO NOT FREEZE. Do not expose the test to direct sunlight.

SPECIMEN COLLECTION AND PREPARATION

The samples from feline feces should be used for this test. The specimens should be tested immediately after collection. If samples are not tested immediately, they should be stored at 2-8°C for 24 hours. For longer storage, freeze at -20°C or below. The right amount of feces to collect is shown in the picture at the top.

PROCEDURE OF THE TEST (LOOK AT THE TOP)

- All reagents must be at room temperature (15-30°C) before use.
- Using 2 swabs, collect the sample from feline feces or directly from the rectum.
- Insert one swab into the specimen tube for FPV/FCoV (P) and the other swab into the specimen tube for Giardia (G). Stir with the swab for at least 10 seconds so the sample and the diluent mix well. Discard the swabs while squeezing the swab against the tube wall.
- Remove the test plate from the foil pouch, and place it on a flat and dry surface.
- Using the disposable pipettes provided, take the supernatant from the specimen tube. Use different pipettes for assay diluent (P) and assay diluent (G). If the feces particles are too big, wait for 1 minute until the large particles are settled, and then take the supernatant of the sample.
- Add four (4) drops into each sample well.
- Read the result after 10 minutes. Do not read after 20 minutes.

INTERPRETATION OF THE TEST (LOOK AT THE TOP)

- Negative result – only the control line (C) appears in the result window
- Positive result – both test line (T) and control line (C) appear in the result window – indicates the presence of target antigens. Even a fade line states a positive result.
- Invalid result – if the control line (C) does not appear, the result might be considered invalid. The sample should be retested.

LIMITATIONS OF THE TEST

Although the Vet Expert FPV/FCoV/Giardia Ag test is very accurate, a low incidence of false results can occur, in case of questionable results, other tests should be performed. Definitive clinical diagnosis should not be based on the results of a single test, but must be an effect of extended clinical and laboratory analysis.

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Sensibilité: FPV 97.1%/FCoV 92.9%/Giardia 99.9%
Spécificité: FPV 98.5%/FCoV 98.1%/Giardia 99.9%

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