

### Product labels and packaging

**1. Product**

PTB Ab ELISA 480

**2. Catalogue number**

PTB480

**3. Outer labels (examples)**

LABEL 1


#### PTB Ab ELISA 480

**REF** PTB480



**LOT** 0100102034

**GMP**

 2025-11



In cooperation with Veterinary Research Institute, Czech Republic

LABEL 2

#### PTB Ab ELISA 480

**MICROPLATE**

5 x 96

**LOT** 0100102037

**NCS**

1 x 0.4 ml

**LOT** 0100102038

**PCS-L**

1 x 0.4 ml

**LOT** 0100102039

**CONJUGATE**

1 x 0.8 ml

**LOT** 0100102040

**DILUENT 14**

1 x 60 ml

**LOT** 0100102041

**CONJ. DILUENT 16**

1 x 60 ml

**LOT** 0100102042

**SUBSTRATE 5**

1 x 60 ml

**LOT** 0100096031

**WASH 20x**

2 x 60 ml

**LOT** 0100099438

**STOP**

1 x 60 ml

**LOT** 0100099552

**LOT** 0100102034

2025-11




**REF** PTB480


V 2




4. Inner labels (examples)

MICROPLATE


 **TestLine** PTB

**MICROPLATE** **LOT** 0100102037


 2025-11


  **GMP** 


NCS

 **TestLine** PTB


**NCS** **LOT** 0100102038

 2025-11


**GMP** 


0.4 ml 


PCS-L

 **TestLine** PTB


**PCS-L** **LOT** 0100102039

 2025-11


**GMP** 


0.4 ml 


CONJUGATE

 **TestLine** PTB

**CONJUGATE** **LOT** 0100102040

 2025-11

**GMP** 

0.8 ml 

DILUENT 14

**B  
G** | **TestLine** PTB

**DILUENT 14**  
60 ml



**LOT** 0100102041

2025-11

**GMP**

2°C 8°C

CONJ. DILUENT 16

**B  
G** | **TestLine** PTB

**CONJ. DILUENT 16**  
60 ml



**LOT** 0100102042

2025-11

**GMP**

2°C 8°C

SUBSTRATE 5

**B  
G** | **TestLine** EIA

**SUBSTRATE 5**  
60 ml



**LOT** 0100096031

2025-12

**GMP**

2°C 8°C

WASH

**B  
G** | **TestLine** EIA

**WASH** 20x

60 ml



**LOT** 0100099438

2026-12

**GMP**

2°C 8°C

STOP

**B  
G** | **TestLine** EIA

☎ +420 549 121 239

**STOP**

60 ml



**LOT** 0100099552

2026-09


**GMP** 2°C 8°C



Contains methanesulphonic acid


TestLine Clinical Diagnostics s.r.o.  
Křižkova, 612 00 Brno-Královo Pole, Czech Republic

**5. Kit box with manufacturer's address**



**TestLine Clinical Diagnostics s.r.o.**  
Křižíkova 68, 612 00 Brno, Czech Republic  
[www.testlinecd.com](http://www.testlinecd.com)

**6. Used symbols**



Temperature limitation



Keep dry



Expiry date



Lot number



Manufactured by



Consult instructions (eIFU)



Catalogue number



Number of tests

**GMP**

Good manufacturing practice